

## Variable Volume Macro Pipette

**MODEL 3100**

- Four different volume models from 1 mL to 10mL.
- Light weight and slender body.
- Easy to operate.
- Liquid touches only disposable tip.
- PTFE piston ensures complete sealing of the syringe.



## Specifications

Cat. No	Volume range(mL)	Increment(mL)
00-31N-1000	0.2 – 1.0	0.1
00-31N-2000	0.4 – 2.0	0.1
00-31N-5000	1.0 – 5.0	0.5
00-31N-10000	2.0 – 10.0	0.5

- Accuracy and Precision  $\pm 1.0\%$  /  $\leq 1.0\%$

## ● Disposable tips

Cat. No	Volume range(mL)	Color	Q'ty(bag)
00-MPT-1	0.2 – 1.0	Clear	1000
00-MPT-2	0.4 – 2.0	Clear	200
00-BMT2-X	1.0 – 5.0	Clear	200
00-BMT2-Z	2.0 – 10.0	Clear	200

## Variable Volume Ultra Macro Pipette

**MODEL 5100**

- Two models, each with three standard volume settings.
- Volume change done by placing a volume spacer in the pipette.
- Any special volume spacer between 15 to 25 mL for 5100S and 30 to 50 mL for 5100L is available for special order.



## Specifications

Cat. No	Volume setting(mL)		
00-5100S	15.0	20.0	25.0
00-5100L	30.0	40.0	50.0

- Accuracy and Precision  $\pm 1.0\%$  /  $\leq 1.0\%$

## ● Disposable tips

Cat. No	Volume range(mL)	Q'ty(bag)
00-51TS	25.0	5
00-51TL	50.0	5