

## Pipette Controller

**Pipette Mate NEO**

- User Friendly pipette controller, 167g only.
- Long-life Lithium Ion battery installed, 4-hours charging for 1700 times.
- Dial-typed Adjustable speed controller for your choice.
- UV resistant for Clean Bench work.
- Self-standing type without holder.
- Fully autoclavable spare parts. ( Nosecone, Sylicon Holder and Filter )
- Complies with CE and RoHS.



## Specifications

Cat. No	Description
00-PMNEO	Pipette Mate NEO
Batteries	Lithium Ion (3.6V, 600mAh)
Charger	100/240V 50/60Hz,0.2A
Filter	20mm×0.2μm (Hydrophobic)
Pipette Types	plastic or glass from 1 to 100mL *Plastics and glass pipettes not sold by Nichiryo

Cat. No	Description
00-PM-N003	Battery cover
00-PM-N005	Head fitting
00-PM-N013	Tube holder
00-PM-N018	Battery L14430-600mAh /w cable
00-PM-N100	Filter, 0.2μm
00-PM-N024	Travel charger 5V US
00-PM-N030	Travel charger base unit
00-PM-N025	Travel charger UK adaptor*
00-PM-N026	Travel charger AU adaptor*
00-PM-N027	Travel charger EU

\*UK and AU adaptors need a travel charger base unit.