

31st May 2023

Nichiryo Co., Ltd.

Dear Customers

Official notification for the suspension of items

We would like to express our sincere gratitude for your continued patronage.

We have confirmed some phenomena that the bottle of the ACCUPENSER JR. which is manufactured and sold by us are broken when it is autoclaved. Therefore, we are afraid that decided to suspend the supplying of ACCUPENSER JR. for 5mL and 10mL until the safety of their quality are confirmed. In addition, we would NOT like you to do autoclave or pour hot liquids into those products that applied following Notification Item. The bottle may crack if you use it in a sudden temperature change. We are currently investigating the cause of this matter. As soon as there is progress, we will promptly disclose information and take appropriate correspondence.

We sincerely apologize for our inconvenience, and your understanding and cooperation would be greatly appreciated.

- Suspension Items : Cat No.: 00-AJ-5
Cat No.: 00-AJ-10
- Suspension Start Date : May 31, 2023

- Notification Item : Cat No.: 00-AJ-5
S/N A18Y02811~A23501261
Date of purchase: November 2018 to May 2023

Cat No.: 00-AJ-10
S/N A22300011~A23501021
Date of purchase: March 2022 to May 2023

- Position of Serial Number:
*See next page

Position of Serial Number

Engraved on top of screw cap



Serial Number