

Pipette Stand



Cat. No 00-NAR-STD
 Storing capacity ... 8pcs
 Compatible models ... Nichipet Air*
 Nichipet Air -Multi Channel-
 Nichipet Premium
 Nichipet Premium LT
 NICHIRYO Le



Cat. No 00-MLT-STD2 (ABS)
 Storing capacity ... 6pcs
 Compatible models ... Nichipet premium LT
 Nichipet Premium
 Nichipet EX II
 Nichipet EX Plus II
 Nichiryo Le
 Nichipet EX II MULTI



Cat. No 00-NLE-STD (ABS)
 Storing capacity ... 1pcs
 Compatible models ... Nichipet premium LT
 Nichipet Premium
 Nichipet EX II
 Nichipet EX Plus II
 Nichiryo Le



Cat. No 00-BM-STD
 (Acrylic resin)
 Storing capacity ... 6pcs
 Compatible models ... Nichipet premium LT
 Nichipet Premium
 Nichipet EX II
 Nichipet EX Plus II
 Nichipet F II & V II

Dispenser Stand



Variable Repetitive Syringe Dispenser MODEL 8100

Cat. No 00-8100S
 (Acrylic resin)
 Storing capacity 2pcs



Ergonomically Designed Repetitive Dispenser Nichimate Stepper

Cat. No 00-NS-STD
 (Acrylic resin)
 Storing capacity 2pcs

Tray



Cat. No 00-88-TY
 Q'ty 50pcs
 Compatible models ... Nichipet EX II MULTI
 MODEL 8800

Tube Opener



Cat. No 00-NTO-001
 Q'ty 3pcs
 Compatible models ... Nichipet premium LT
 Nichipet Premium
 Nichipet EX II
 Nichipet EX Plus II

Volume adjustment tool



Cat. No 00-NAR-AT
 Q'ty 1pcs
 Compatible models ... Nichipet Air
 Nichipet Air -Multi Channel-

Filter

For NPP-1000
 Cat. No 00-NP-0601000
 Q'ty 10pcs



For NPLO II -1000
 Cat. No 00-LO2-0601000
 Q'ty 10pcs



For NPX II -1000
 Cat. No 00-NX2-0601000
 Q'ty 10pcs



For NPP,NPX II,NPLO II -5000,10000
 Cat. No ... 00-NX2-0615000
 Q'ty 100pcs



For NAR-5000,10000
 Cat. No ... 00-NAR-F15000
 Q'ty 10pcs



CUSTOM MADE & OEM

Custom-Made products

Custom-made Pipettes

- Talk to us when you wish to have pipettes or molding products that you cannot find in the market.
- Full Turnkey Contract starting from planning, developing, manufacturing and maintaining products will be possible.
- Modification to our existing product is also welcomed.
- Our custom-made molding products include disposable tip, micro plate, cell culture vessels and other plastic product.



Custom-Made Automatic Dispenser

NSS-1000 II Basic

- Cost effective custom made automatic dispensing systems
- Stand alone dispensers as well as dispensing components to be integrated into larger system units will be possible.
- Working areas ranging from A4-sized desk top type to larger work station type will be possible.
- Software to meet customers' requirements concerning C-language and Visual Basic will be possible.



Custom-made system example

Appearance	Name	Main features	Appearance	Name	Main features
	Dilution system for mechanical oil	Function to prevent contamination with air pressure feed-in. 1000mm(W)×500mm(D)×500mm		Sample preparation system by disposable tip	Simple dispensing of samples from test tubes to micro plates using disposable tips. 500mm(W)×800mm(D)×150mm
	Two channel micro plate dilution system	Dilutes samples on micro plates. Two channels with disposable and washable tips. 700mm(W)×450mm(D)×700mm		Automatic picking system for Electrophoresis gel	Determines sample spots fractionated by electrophoresis through image analysis and performs safe picking with the rotation and suction feature. 756mm(W)×700mm(D)×750mm
	Multiple dispensing and dilution system	Automatic dispenser with two syringe pumps for dispensing processes and liquid line cleaning. 800mm(W)×450mm(D)×700mm		Five channel high-speed dispensing system	System for dispensing five different liquids using five syringe pumps and five probes. 800mm(W)×500mm(D)×700mm
	Dispenser for PCR preparation	Reformatting from four 96-well micro plates to one 384-well micro plate. 500mm(W)×400mm(D)×500mm		Powder sample dispensing system	High accuracy dispenses of small powder samples with the suction and discharging function. Can also be connected to electronic balances or combined with liquid dispensing. 720mm(W)×700mm(D)×550mm
	Five channel high throughput dilution system	Automatic tip interval span feature and stacker for tips and micro plates are included. 1000mm(W)×610mm(D)×500mm			