

Disposable Tips

Premium Tip

Tip

- The fine design and finishing of the top of the tips reduce the surface tension of liquid and ensure the tip to be residue free.
- The excellent sealing and discharge ability enable high accuracy of pipetting.
- Chemical and heat resistant, for which the advanced material "Highly purified polypropylene" is used.
- Fully autoclavable. (121°C for 20 minutes) Racks can be autoclaved as well. (*Exceeding the recommended standard for autoclaving may cause deformities)
- The tips being able to be attached directly from the racks avoids dirt and dust from getting in, thus rising the efficiency of the operations.
- DNase/RNase free (excludes X/Z series)
- UV resistant




Specifications

* The tips are available for pipette models with in the volume range of the tip mentioned above.

Type	Cat. No.	Specification	Volume range(μL)	Color	Q'ty(box)	Usable models
UT	00-BMT2-UT	Bulk tip	0.1~10	Natural	1000	NAR-2, 10 / NAR-10M8, 10M12 NPP-2 / NLT-2 NPX-NPX2-2 / NPLO-NPLO2-2 NP7-8V, 12V / NPM-8V, 12V
	00-BMT2-UTWR	Racked tip			960 (96pcs.×10 ² /s)	
	00-BMT2-UTWRB	Racked tip (sterilized)			960 (96pcs.×10 ² /s)	
SSW	00-BMT2-SSW	Bulk tip	0.5~10	Natural	1000	NPP-10 / NLT-10 NPX-NPX2-10 / NPLO-NPLO2-10 NLE-10 / NPF2-10
	00-BMT2-SSWR	Racked tip			960 (96pcs.×10 ² /s)	
	00-BMT2-SSWRB	Racked tip (sterilized)			960 (96pcs.×10 ² /s)	
SG	00-BMT2-SG	Bulk tip	2~200	Natural	1000	NAR-20, 100, 200 / NAR-100M8, 100M12, 200M8, 200M12 NPP-20, 100, 200 / NLT-20, 100, 200 NPX-20, 100, 200 / NPLO2-20, 100, 200 NLE-20, 100, 200 / NPF2-20, 100, 200 NPV2-S, M / NP7-8S, 8L, 12S, 12L NPM-8S, 8L, 12S, 12L
	00-BMT2-SGR	Racked tip			960 (96pcs.×10 ² /s)	
	00-BMT2-SGRB	Racked tip (sterilized)			960 (96pcs.×10 ² /s)	
SE	00-BMT2-SE	Bulk tip	2~200	Natural	1000	NPX-20, 100, 200 NPLO-20, 100, 200
	00-BMT2-SER	Racked tip			960 (96pcs.×10 ² /s)	
	00-BMT2-SERB	Racked tip (sterilized)			960 (96pcs.×10 ² /s)	
K	00-BMT2-K	Bulk tip	30~300	Natural	1000	NAR-300M8, 300M12 NP7-8K, 12K NPM-8K, 12K
	00-BMT2-KR	Racked tip			960 (96pcs.×10 ² /s)	
	00-BMT2-KRB	Racked tip (sterilized)			960 (96pcs.×10 ² /s)	
LG	00-BMT2-LG	Bulk tip	100~1000	Natural	1000	NAR-1000 NPP-1000 / NLT-1000 NPX2-1000 / NPLO2-1000 NLE-1000 / NPF2-500, 1000 / NPV2-L
	00-BMT2-LGR	Racked tip			960 (96pcs.×10 ² /s)	
	00-BMT2-LGRB	Racked tip (sterilized)			960 (96pcs.×10 ² /s)	
LE	00-BMT2-LE	Bulk tip	100~1000	Natural	1000	NPX-1000 NPLO-1000
	00-BMT2-LER	Racked tip			960 (96pcs.×10 ² /s)	
	00-BMT2-LERB	Racked tip (sterilized)			960 (96pcs.×10 ² /s)	
X	00-BMT2-X	Bulk tip	1000~5000	Natural	200	NAR-5000 NPP-5000 NPX-NPX2-5000 NPLO-NPLO2-5000/11DG-5000
	00-BMT2-XR	Racked tip			100 (50pcs.×2 ² /s)	
	00-BMT2-XRB	Racked tip (sterilized)			100 (50pcs.×2 ² /s)	
Z	00-BMT2-Z	Bulk tip	1000~10000	Natural	200	NPP-10000 NPX-NPX2-10000 NPLO-NPLO2-10000
	00-BMT2-ZR	Racked tip			80 (40pcs.×2 ² /s)	
	00-BMT2-ZRB	Racked tip (sterilized)			80 (40pcs.×2 ² /s)	
ZS	00-BMT2-ZS	Bulk tip	1000~10000	Natural	200	NAR-10000
	00-BMT2-ZSR	Racked tip			80 (40pcs.×2 ² /s)	
	00-BMT2-ZSRB	Racked tip (sterilized)			80 (40pcs.×2 ² /s)	

TIP

(measure : mm)



Filter Tip

- Prevents cross-contamination in PCR testing.
- Suitable for virus/bacteria and RI research.
- Sterilized racked tips available for DNA free research.
- No effect on the accuracy by the inside filter.
- Efficient for avoiding carry over by cellulose acetate filters.
- DNase/RNase free



Specifications

Type	Cat. No.	Specification	Volume range(μL)	Color	Q'ty(box)	Usable models
UT	00-BMT2-FUT	Bulk tip	0.1~10	Natural	1000	NAR-2, 10 / NAR-10M8, 10M12 NPP-2 / NLT-2 NPX-NPX2-2 / NPLO-NPLO2-2 NP7-8V, 12V / NPM-8V, 12V
	00-BMT2-FUTRB	Racked tip (sterilized)			960 (96pcs.×10 ² /s)	
SS	00-BMT2-FSS	Bulk tip	0.5~10	Natural	1000	NAR-10 NPP-10 / NLT-10 NPX-NPX2-10 / NPLO-NPLO2-10 NLE-10 / NPF2-10
	00-BMT2-FSSRB	Racked tip (sterilized)			960 (96pcs.×10 ² /s)	
SG	00-BMT2-FSG	Bulk tip	2~20	Natural	1000	NAR-20 NPP-20 / NLT-20 NPX-NPX2-20 / NPLO-NPLO2-20 NLE-20 / NPF2-20
	00-BMT2-FSGRB	Racked tip (sterilized)			960 (96pcs.×10 ² /s)	
SE	00-BMT2-FSE	Bulk tip	10~100	Natural	1000	NAR-100 / NAR-100M8, 100M12 NPP-100 / NLT-100 NPX-NPX2-100 / NPLO-NPLO2-100 NLE-100 / NPF2-100
	00-BMT2-FSERB	Racked tip (sterilized)			960 (96pcs.×10 ² /s)	
AG	00-BMT2-FAG	Bulk tip	20~200	Natural	1000	NAR-200 / NAR-200M8, 200M12 NPP-200 / NLT-200 / NPX-NPX2-200 NPLO-NPLO2-200 / NPF2-200 NPV-M / NPM-8LP, 12LP
	00-BMT2-FAGRB	Racked tip (sterilized)			960 (96pcs.×10 ² /s)	
LG	00-BMT2-FLG	Bulk tip	100~1000	Natural	1000	NAR-1000 NPP-1000 / NLT-1000 NPX-NPX2-1000 / NPLO-NPLO2-1000 NLE-1000 / NPF2-1000 / NPV2-L
	00-BMT2-FLGRB	Racked tip (sterilized)			960 (96pcs.×10 ² /s)	
LE	00-BMT2-FLE	Bulk tip	100~1000	Natural	1000	NPX-NPX2-1000 NPLO-NPLO2-1000 / NLE-1000 NPF2-1000 / NPV2-L
	00-BMT2-FLERB	Racked tip (sterilized)			960 (96pcs.×10 ² /s)	

Smart Pack

- Easy & quick tip setting
- Eco-friendly & low-waste

How to refill

1. Lift up the plate of the top deck of tips.
2. Holding the empty rack in one hand and the plate in the other, align the first row of tips in the rack. Insert the remaining tips into the rack.
3. Release the plate.



Specifications

Type	Cat. No.	Specification	Volume range(μL)	Color	Q'ty(box)	Usable models
Smart Pack UT	00-ETS-JUTW	Refill tip	0.1~10	Natural	96pcs×10decks	NPP-2 / NLT-2 NPX-NPX2-2 / NPLO-NPLO2-2 NP7-8V, 12V / NPM-8V, 12V
	00-ETS-JUTW10	Rack			96pcs×10decks×10boxes	
	00-ETS-UTWRK	Rack			5racks	
Smart Pack SG	00-ETS-JSG	Refill tip	2~200	Natural	96pcs×5decks×2boxes	NPP-20, 100, 200 / NLT-20, 100, 200 / NPX2-20, 100, 200 NPLO2-20, 100, 200 / NLE-20, 100, 200 NPF2-20~200 / NPV2-S, M / NP7-8S, 8L, 12S, 12L NPM-8SP, 8LP, 12SP, 12LP
	00-ETS-JSG10	Rack			96pcs×5decks×2×10boxes	
	00-ETS-SGRK	Rack			5racks	
Smart Pack LG	00-ETS-JLG	Refill tip	100~1000	Natural	96pcs×5decks×2boxes	NPP-1000 / NLT-1000 NPX2-1000 / NPLO2-1000 NLE-1000 / NPF2-500, 1000 / NPV2-L
	00-ETS-JLG10	Rack			96pcs×5decks×2×10boxes	
	00-ETS-LGRK	Rack			5racks	