

Standard Tip

NichimateAQUATip

Bulk tip , Refill tip

■ The N10-RF and N200-RF are compatible with the following racks: N10RK, N10LRK, and N200RK.



Refill tip



Bulk tip



Type	Length	Cat No.	Specification	Q'ty(box)
10 μ L	31.3mm	00-AQA-N10	Bulk tip	1,000pcs
		00-AQA-N10-RF	Refill tip	960pcs (96pcs \times 10stacks)
10 μ L (Long)	48.0mm	00-AQA-N10L	Bulk tip	1,000pcs
200 μ L	53.0mm	00-AQA-N200	Bulk tip	1,000pcs
		00-AQA-N200-RF	Refill tip	960pcs (96pcs \times 5stacks \times 2)
1000 μ L	81.9mm	00-AQA-N1000	Bulk tip	1,000pcs
		00-AQA-N1000-RF	Refill tip	960pcs (96pcs \times 5stacks \times 2)

Rack (Made of Transparent polycarbonate)

■ This rack is exclusively designed for Nichimate Aqua Tips and is made of polycarbonate. Tips are not included.



Type	Cat No.	Q'ty(box)
10 μ L	00-AQA-N10RK (Orange●)	5 pcs.
10 μ L (Long)	00-AQA-N10LRK (Red●)	
200 μ L	00-AQA-N200RK (Yellow●)	
1000 μ L	00-AQA-N1000RK (Blue●)	



Standard Tip (Sterilized Raked Tip)

NichimateAqua MarineTip



10 μ L

10 μ L (Long)

200 μ L

200 μ L (Long)

1000 μ L

Type	Length	Cat No.	Specification	Q'ty(box)
10 μ L	32mm	00-MRN-10RB	Racked tip	960pcs (96pcs \times 10 $^{\circ}$ /s)
10 μ L (Long)	46mm	00-MRN-10LRB		
200 μ L	51mm	00-MRN-200RB		
200 μ L (Long)	90mm	00-MRN-200LRB		
1000 μ L	88mm	00-MRN-1000RB		

● This product is not recommended for repeated use with autoclave sterilization.

Tips with Excellent Liquid Drainage

Low Retention Tip



- Standard tips with high transparency and excellent liquid drainage.
- Especially effective for handling viscous liquids and surfactants that generate foam during dispensing.
- Ideal for preparing PCR reagents due to the above features.

- Bulk and Refill are non-sterile; Rack tips are provided sterile.
- DNase/RNase/Pyrogen Free.
- For research use only. Not for use in medical or clinical applications.

Specifications

Specification	Type	Length	Cat No.	Q'ty(box)
Bulk tip	10 μ L	31.1mm	00-LRT-10	1,000pcs
	10 μ L (Long)	45.8mm	00-LRT-10L	
	200 μ L	50.8mm	00-LRT-200	
	300 μ L	55.9mm	00-LRT-300	
	1250 μ L	90.2mm	00-LRT-1250	
Refill tip	10 μ L	31.1mm	00-LRT-10-RF	960pcs (96pcs \times 5stacks \times 2)
	10 μ L (Long)	45.8mm	00-LRT-10L-RF	
	200 μ L	50.8mm	00-LRT-200-RF	
	300 μ L	55.9mm	00-LRT-300-RF	768pcs (96pcs \times 4stacks \times 2)
	1250 μ L	90.2mm	00-LRT-1250-RF	
Racked tip	10 μ L	31.1mm	00-LRT-10RB	960pcs (96pcs \times 10 $^{\circ}$ /s)
	10 μ L (Long)	45.8mm	00-LRT-10LRB	
	200 μ L	50.8mm	00-LRT-200RB	
	200 μ L (Long)	91.1mm	00-LRT-200LRB	192pcs (96pcs \times 2 $^{\circ}$ /s)
	300 μ L	55.9mm	00-LRT-300RB	960pcs (96pcs \times 10 $^{\circ}$ /s)
	1250 μ L	90.2mm	00-LRT-1250RB	

● The racks can be autoclaved (121°C, 20 minutes); however, but should be replaced if deformation or discoloration occurs (Nichiryo does not guarantee that deformation or discoloration will not occur). ● Bulk and refill packages cannot be autoclaved as is. Please transfer the tips to a compatible rack, glass container, or similar ones before autoclaving. ● Autoclaving should be performed according to your laboratory manual or under the guidance of an experienced person.

Filtered Tips with Excellent Liquid Drainage

Low Retention Filtered Tip

- Features a hydrophobic HDPE filter to minimize aerosol contamination inside the micropipette during aspiration.
- DNase/RNase/Pyrogen Free.

- For research use only. Not for use in medical or clinical applications.

Specifications

Type	Length	Cat No.	Specification	Q'ty(box)
10 μ L	31.1mm	00-LRT-F10RB	Racked tip (sterilized)	960pcs (96pcs \times 10 $^{\circ}$ /s)
10 μ L (Long)	45.8mm	00-LRT-F10LRB		
20 μ L	50.8mm	00-LRT-F20RB		
100 μ L		00-LRT-F100RB		192pcs (96pcs \times 2 $^{\circ}$ /s)
200 μ L		00-LRT-F200RB		
200 μ L (Long)	91.1mm	00-LRT-F200LRB		
300 μ L	55.9mm	00-LRT-F300RB	960pcs (96pcs \times 10 $^{\circ}$ /s)	
1250 μ L	90.2mm	00-LRT-F1250RB		

● The racks can be autoclaved (121°C, 20 minutes), but should be replaced if warped or discolored.