

Standard Tip (Racked, Sterilized)

## NichimateAqua *Marine Tip*

DNase  
RNase  
**FREE**

Pyrogen  
**FREE**



- Reasonably priced and has good liquid drainage performance
- 96 tips per rack, supplied in 10 racks
- Electron beam sterilized, DNase/RNase/Pyrogen (endotoxin) free



10µL



10µL (Long)




200µL



200µL (Long)



1000µL



### SPECIFICATIONS

Package	Volume	Cat. No.	Length of Tip	Q'ty (Box)	Sterilization
Racked	10µL	00-MRN-10RB	32mm	96 pcs. × 10 c/s	E. beam
	10µL (Long)	00-MRN-10LRB	46mm		
	200µL	00-MRN-200RB	51mm		
	200µL (Long)	00-MRN-200LRB	90mm		
	1000µL	00-MRN-1000RB	88mm		

● Please ask us the compatibility with the micropipette models in advance. ● We do not recommend repeated use of this product with autoclaving. ● The packing box of this product is for protection during transportation, and may have some scratches or dirt. Please be noted that we cannot replace the packing box or product if the product inside is not affected. Your understanding and cooperation would be appreciated to reduce unnecessary consumption of resources. ● Please be noted that the product standards and specifications are subject to change without notice for improvement purpose.

# NICHIRYO CO.,LTD.

<https://www.nichiryo.co.jp/en>

[info@nichiryo.co.jp](mailto:info@nichiryo.co.jp)

● For inquiries, please use the inquiry form on our website or email us at the above address. ● Published on April 1, 2025.