


Pipette Controller

Pipette Mate NEO

- User Friendly pipette controller, 167g only
- Long-life Lithium Ion battery is installed. About 1700 times dispensing are possible by full charge (4 hours).
- Dial-typed adjustable speed controller for your choice.
- UV resistant for clean bench work.
- Self- standing type without holder.
- Fully autoclavable spare parts. (Nosecone, Sylicon Holder, and Filter)
- Complies with CE and RoHS.



【Specifications】



Specifications

Cat. No	Description
00-PMNEO	Pipette Mate NEO
Batteries	Lithium Ion (3.6V, 600mAh)
Charger	100/240V 50/60Hz, 0.2A
Filter	20mm×0.2μm (Hydrophobic)
Pipette Types	Plastic or glass from 1 to 100mL * Plastics and glass pipettes not sold by Nichiryo

Cat. No	Description
00-PM-N005	Head fitting
00-PM-N013	Tube holder
00-PM-N100	Filter, 0.2μm
00-PM-N024	Travel charger 5V US
00-PM-N030	Travel charger base unit
00-PM-N025	Travel charger UK adaptor*
00-PM-N026	Travel charger AU adaptor*
00-PM-N027	Travel charger EU

* UK and AU adaptors need a travel charger base unit.