

Disposable Tips

Premium Tip

[Except X/Z/ZS]





- The fine design and finishing of the top of the tips reduce the surface tension of liquid and ensure the tip to be residue free.
- The excellent sealing and discharge ability enable high accuracy of pipetting.
- Chemical and heat resistant, for which the advanced material "Highly purified polypropylene" is used.
- Fully autoclavable. (121°C for 20 minutes) Racks can be autoclaved as well. (*Exceeding the recommended standard for autoclaving may cause deformities)
- The tips being able to be attached directly from the racks avoids dirt and dust from getting in, thus rising the efficiency of the operations.
- DNase/RNase free (excludes X/Z series)
- UV resistant

Specifications

*Usable models refer to the micropipettes calibrated by our company using this tips. For compatibility with micropipettes from other manufacturers, please contact us.

*The tips are available for pipette models with in the volume range of the tip mentioned above.

Type		Cat. No.	Specification	Volume range(μL)	Color	Q'ty(box)	Usable models
UT		00-BMT2-UT	Bulk tip	0.1~10	Natural	1,000pcs	NAR-2, 10 / NAR-10M8, 10M12
		00-BMT2-UTWR	Racked tip			960pcs (96pcs×10 ² /s)	NPX3-2, 10 / NPP-2 / NLT-2
		00-BMT2-UTWRB	Racked tip (sterilized)				NPX2-2 / NPLO-NPLO2-2
SSW		00-BMT2-SSW	Bulk tip	0.5~10		1,000pcs	NPX3-2, 10
		00-BMT2-SSWR	Racked tip			960pcs (96pcs×10 ² /s)	NPP-10 / NLT-10 / NPX3-10
		00-BMT2-SSWRB	Racked tip (sterilized)				NPX2-10 / NPLO-NPLO2-10
SG		00-BMT2-SG	Bulk tip	2~200		1,000pcs	NAR-20, 100, 200 / NAR-100M8, 100M12, 200M8, 200M12
		00-BMT2-SGR	Racked tip			960pcs (96pcs×10 ² /s)	NPX3-20, 100, 200 / NPP-20, 100, 200 / NLT-20, 100, 200
		00-BMT2-SGRB	Racked tip (sterilized)				NPX2-20, 100, 200 / NPLO2-20, 100, 200
SE		00-BMT2-SE	Bulk tip	2~200		1,000pcs	NLE-20, 100, 200 / NPF2-20~200
		00-BMT2-SER	Racked tip			960pcs (96pcs×10 ² /s)	NPV2-S, M / NPM-SS, SL, 12S, 12L
		00-BMT2-SERB	Racked tip (sterilized)				
K		00-BMT2-K	Bulk tip	30~300		1,000pcs	
		00-BMT2-KR	Racked tip			960pcs (96pcs×10 ² /s)	NPX-20, 100, 200
		00-BMT2-KRB	Racked tip (sterilized)				NPLO-20, 100, 200
LG		00-BMT2-LG	Bulk tip	100~1000		1,000pcs	NAR-300M8, 300M12
		00-BMT2-LGR	Racked tip			960pcs (96pcs×10 ² /s)	NPM-8K, 12K
		00-BMT2-LGRB	Racked tip (sterilized)				
LE		00-BMT2-LE	Bulk tip	100~1000		1,000pcs	NAR-1000 / NPX3-1000
		00-BMT2-LER	Racked tip			960pcs (96pcs×10 ² /s)	NPP-1000 / NLT-1000
		00-BMT2-LERB	Racked tip (sterilized)				NPX2-1000 / NPLO2-1000
X		00-BMT2-X	Bulk tip	1000~5000		200pcs	NLE-1000 / NPF2-500, 1000 / NPV2-L
		00-BMT2-XR	Racked tip			100pcs (50pcs×2 ² /s)	NPX-1000
		00-BMT2-XRB	Racked tip (sterilized)				NPLO-1000
Z		00-BMT2-Z	Bulk tip	1000~10000		200pcs	NAR-5000 / NPX3-5000
		00-BMT2-ZR	Racked tip			80pcs (40pcs×2 ² /s)	NPP-5000
		00-BMT2-ZRB	Racked tip (sterilized)				NPX2-5000
ZS		00-BMT2-ZS	Bulk tip	1000~10000		200pcs	NPLO-NPLO2-5000
		00-BMT2-ZSR	Racked tip			80pcs (40pcs×2 ² /s)	NPP-10000
		00-BMT2-ZSRB	Racked tip (sterilized)				NPX2-10000

TIP

(measure:mm)

