

Autoclavable Bottle-top Dispenser

DISPET EX

- Autoclavable. (121°C for 20 minutes)
- 3 sized adapters for fitting with various types of bottles.
- Chemical-guard materials used for liquid contact parts.
- Adjustable Inlet tube for various types of bottles.
- Liquid returning mechanism for minimizing the waste of reagents.

*Bottles are available as optional parts.

121°C
20min
Autoclavable
γ-ray irradiation
50kGy



Specifications

Cat. No	Volume range (mL)	Increment (mL)	Accuracy (%)	Precision (%)	Thread	Attached adapters
00-DPX-25	0.25 – 2.5	0.05	±0.6	≤0.1	GL-32	GL-28, GL-38, GL-45
00-DPX-50	0.5 – 5.0	0.1	±0.5	≤0.1	GL-32	GL-28, GL-38, GL-45
00-DPX-100	1.0 – 10.0	0.2	±0.5	≤0.1	GL-32	GL-28, GL-38, GL-45
00-DPX-250	2.5 – 25.0	0.5	±0.5	≤0.1	GL-45	GL-28, GL-32, GL-38
00-DPX-500	5.0 – 50.0	1	±0.5	≤0.1	GL-45	GL-28, GL-32, GL-38
00-DPX-1000	10.0 – 100.0	2	±0.5	≤0.1	GL-45	GL-28, GL-32, GL-38

Adapter

Cat. No	Size	Cat. No	Size	Cat. No	Size
00-DPX-003A	28mm	00-DPX-003F	38mm	00-DPX-003D	40mm
00-DPX-003B	32mm	00-DPX-003C	38mm	00-DPX-003E	45mm

00-DPX-003F...For DPX-25,50,100
00-DPX-003C...For DPX-250,500,1000

Bottle-top Dispenser

DISPET Mini

- Suitable for repetitively micro volume dispensing.
- Light and compact model for easy-to-carry-out.
- Chemical-guard materials used for liquid contact parts.

*Bottles are available as optional parts.



Specifications

Cat. No	Volume range (μL)	Accuracy (%)	Precision (%)	Thread
00-DPM-250	100 / 250	±2.0	≤0.4	GL-28
00-DPM-1000	500 / 1000	±1.0	≤0.2	GL-28
00-DPM-2000	1000 / 2000	±1.0	≤0.2	GL-28