

Hand-held Dispenser

# MODEL 2100DG

- Dispensing from most types of liquid containers is possible.
- Easy to use with fixed volume adjustment.
- Light weight design for one-hand operation.
- The glass cylinder and plunger are protected by plastic housings.



## Specifications

Cat. No	Volume range(mL)	Increment(mL)	Volume readout
00-21DG-1	0.2 – 1.0	0.01	0/20 – 1/00
00-21DG-2	0.4 – 2.0	0.01	0/40 – 2/00
00-21DG-5	1.0 – 5.0	0.02	1/00 – 5/00
00-21DG-10	2.0 – 10.0	0.1	02/0 – 10/0

\* Liquid contact parts ; PTFE, Tygon® Polypropylene, Pyrex® Glass, Ceramic and Platinum Iridium.

## Accuracy and Precision

±1.0% / ≤0.2%