

PIPETTE

Disposable Tips




Premium Tip

 Tip

- The fine design and finishing of the top of the tips reduce the surface tension of liquid and ensure the tip to be residue free.
- The excellent sealing and discharge ability enable high accuracy of pipetting.
- Chemical and heat resistant, for which the advanced material "Highly purified polypropylene" is used.
- Fully autoclavable. (121°C for 20 minutes) Racks can be autoclaved as well. (*Exceeding the recommended standard for autoclaving may cause deformities)
- The tips being able to be attached directly from the racks avoids dirt and dust from getting in, thus rising the efficiency of the operations.
- DNase/RNase free. (excludes X/Z series)
- UV resistant.



Specifications

Type		Cat. No.	Specification	Volume range(μL)	Color	Q'ty(box)	Usable models
UT		00-BMT2-UT	Bulk tip	0.1~10	Natural	1000	NPP-2 / NLT-2 NPX·NPX2-2 / NPLO·NPLO2-2 NP7-8V, 12V / NPM-8VP, 12VP
		00-BMT2-UTWR	Racked tip			960 (96pcs.x10 ⁰ /s)	
		00-BMT2-UTWRB	Racked tip(sterilized)				
SS		00-BMT2-SS	Bulk tip	0.5~10	Natural	1000	NPP-10 / NLT-10 NPX·NPX2-10 / NPLO·NPLO2-10 NLE-10 / NPF2-10
		00-BMT2-SSR	Racked tip			960 (96pcs.x10 ⁰ /s)	
		00-BMT2-SSRB	Racked tip(sterilized)				
SG		00-BMT2-SG	Bulk tip	2~200	Natural	1000	NPP-20, 100, 200 / NLT-20, 100, 200 NPX2-20, 100, 200 / NPLO2-20, 100, 200 NLE-20, 100, 200 / NPF2-20~200 NPV2-S, M / NP7-8S, 8L, 12S, 12L NPM-8SP, 8LP, 12SP, 12LP
		00-BMT2-SGR	Racked tip			960 (96pcs.x10 ⁰ /s)	
		00-BMT2-SGRB	Racked tip(sterilized)				
SE		00-BMT2-SE	Bulk tip	2~200	Natural	1000	NPX-20, 100, 200 NPLO-20, 100, 200
		00-BMT2-SER	Racked tip			960 (96pcs.x10 ⁰ /s)	
		00-BMT2-SERB	Racked tip(sterilized)				
K		00-BMT2-K	Bulk tip	30~300	Natural	1000	NP7-8K, 12K NPM-8KP, 12KP
		00-BMT2-KR	Racked tip			960 (96pcs.x10 ⁰ /s)	
		00-BMT2-KRB	Racked tip(sterilized)				
LG		00-BMT2-LG	Bulk tip	100~1000	Natural	1000	NPP-1000 / NLT-1000 NPX2-1000 / NPLO2-1000 NLE-1000 / NPF2-500, 1000 / NPV2-
		00-BMT2-LGR	Racked tip			960 (96pcs.x10 ⁰ /s)	
		00-BMT2-LGRB	Racked tip(sterilized)				
LE		00-BMT2-LE	Bulk tip	100~1000	Natural	1000	NPX-1000 NPLO-1000
		00-BMT2-LER	Racked tip			960 (96pcs.x10 ⁰ /s)	
		00-BMT2-LERB	Racked tip(sterilized)				
X		00-BMT2-X	Bulk tip	1000~5000	Natural	200	NPP-5000 NPX·NPX2-5000 NPLO·NPLO2-5000 11DG-5000
		00-BMT2-XR	Racked tip			100 (50pcs.x2 ⁰ /s)	
		00-BMT2-XRB	Racked tip(sterilized)				
Z		00-BMT2-Z	Bulk tip	1000~10000	Natural	200	NPP-10000 NPX·NPX2-10000 NPLO·NPLO2-10000
		00-BMT2-ZR	Racked tip			80 (40pcs.x2 ⁰ /s)	
		00-BMT2-ZRB	Racked tip(sterilized)				

*The tips are available for pipette models with in the volume range of the tip mentioned above.

TIP

