

PIPETTE

Disposable Tips

Premium Tip

Tip

- The fine design and finishing of the top of the tips reduce the surface tension of liquid and ensure the tip to be residue free.
- The excellent sealing and discharge ability enable high accuracy of pipetting.
- Chemical and heat resistant, for which the advanced material "Highly purified polypropylene" is used.
- Fully autoclavable. (121°C for 20 minutes) Racks can be autoclaved as well. (*Exceeding the recommended standard for autoclaving may cause deformities)
- The tips being able to be attached directly from the racks avoids dirt and dust from getting in, thus rising the efficiency of the operations.
- DNase/RNase free. (excludes X/Z series)
- UV resistant.



Specifications

* The tips are available for pipette models with in the volume range of the tip mentioned above.

Type	Cat. No.	Specification	Volume range(μL)	Color	Q'ty(box)	Usable models
UT	00-BMT2-UT	Bulk tip	0.1~10	Natural	1000	NAR-2, 10 NPP-2 / NLT-2 NPX·NPX2-2 / NPLO·NPLO2-2 NP7-8V, 12V / NPM-8V, 12V
	00-BMT2-UTWR	Racked tip			960	
	00-BMT2-UTWRB	Racked tip (sterilized)			(96pcs. × 10 ² /s)	
SSW	00-BMT2-SSW	Bulk tip	0.5~10	Natural	1000	NAR-10 NPP-10 / NLT-10 NPX·NPX2-10 / NPLO·NPLO2-10 NLE-10 / NPF2-10
	00-BMT2-SSWR	Racked tip			960	
	00-BMT2-SSWRB	Racked tip (sterilized)			(96pcs. × 10 ² /s)	
SG	00-BMT2-SG	Bulk tip	2~200	Natural	1000	NAR-20, 100, 200 NPP-20, 100, 200 / NLT-20, 100, 200 NPX2-20, 100, 200 / NPLO2-20, 100, 200 NLE-20, 100, 200 / NPF2-20~200 NPV2-S, M / NP7-8S, 8L, 12S, 12L NPM-8S, 8L, 12S, 12L
	00-BMT2-SGR	Racked tip			960	
	00-BMT2-SGRB	Racked tip (sterilized)			(96pcs. × 10 ² /s)	
SE	00-BMT2-SE	Bulk tip	2~200	Natural	1000	NPX-20, 100, 200 NPLO-20, 100, 200
	00-BMT2-SER	Racked tip			960	
	00-BMT2-SERB	Racked tip (sterilized)			(96pcs. × 10 ² /s)	
K	00-BMT2-K	Bulk tip	30~300	Natural	1000	NP7-8K, 12K NPM-8K, 12K
	00-BMT2-KR	Racked tip			960	
	00-BMT2-KRB	Racked tip (sterilized)			(96pcs. × 10 ² /s)	
LG	00-BMT2-LG	Bulk tip	100~1000	Natural	1000	NAR-1000 NPP-1000 / NLT-1000 NPX2-1000 / NPLO2-1000 NLE-1000 / NPF2-500, 1000 / NPV2-L
	00-BMT2-LGR	Racked tip			960	
	00-BMT2-LGRB	Racked tip (sterilized)			(96pcs. × 10 ² /s)	
LE	00-BMT2-LE	Bulk tip	100~1000	Natural	1000	NPX-1000 NPLO-1000
	00-BMT2-LER	Racked tip			960	
	00-BMT2-LERB	Racked tip (sterilized)			(96pcs. × 10 ² /s)	
X	00-BMT2-X	Bulk tip	1000~5000	Natural	200	NAR-5000 NPP-5000 NPX·NPX2-5000 NPLO·NPLO2-5000/11DG-5000
	00-BMT2-XR	Racked tip			100	
	00-BMT2-XRB	Racked tip (sterilized)			(50pcs. × 2 ² /s)	
Z	00-BMT2-Z	Bulk tip	1000~10000	Natural	200	NAR-10000 NPP-10000 NPX·NPX2-10000 NPLO·NPLO2-10000
	00-BMT2-ZR	Racked tip			80	
	00-BMT2-ZRB	Racked tip (sterilized)			(40pcs. × 2 ² /s)	
ZS	00-BMT2-ZS	Bulk tip	1000~10000	Natural	200	NAR-10000
	00-BMT2-ZSR	Racked tip			80	
	00-BMT2-ZSRB	Racked tip (sterilized)			(40pcs. × 2 ² /s)	

TIP

(measure : mm)

